Instruction Manual Advanced Series

Stirrers

Hotplates

Hotplate-Stirrers

TABLE OF CONTENTS

Package Contents	1
Warranty	1
Installation	2
Maintenance & Servicing	2
Intended Use	2
Equipment Disposal	2
Environmental Conditions	2
Safety Instructions	3
Standards & Regulations	3
10 x 10 Specifications	4
18 x 18 Specifications	5
25 x 25 Specifications	6
Heating Operating Instructions	7-9
Stirring Operating Instructions	10
Technical Service	11
Troubleshooting	11
Replacement Parts	12-14

PACKAGE CONTENTS

- Hotplate, Stirrer or Hotplate-Stirrer
- 234cm detachable power cord
- Spin bar (stirring models only)
- Instruction manual

WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

INSTALLATION

Upon receiving the VWR Hotplate/Stirrer/Hotplate-Stirrer, check to ensure that no damage has occurred in shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the Hotplate/Stirrer/Hotplate-Stirrer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of 15,2 centimeters from vertical surfaces. Always place the unit on a sturdy work surface.

The Hotplate/Stirrer/Hotplate-Stirrer is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The unit plugs into a 230 volt, 50/60 Hz source.

MAINTENANCE & SERVICING

The Hotplate/Stirrer/Hotplate-Stirrer is built for long, trouble-free, dependable service. It needs no user maintenance beyond keeping the surfaces clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly after the unit has cooled down. Do not use a cleaning agent or solvent on the front panel or top plate which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

CLEANING CERAMIC TOPS:

First remove any burnt-on deposits or spills from the top plate with a scraper (similar to scraping paint off of windowpanes in your home). For your safety, please wear an insulated mit when using a metal scraper. When the top plate has cooled, apply a few dabs of a non-abrasive cleaner over the surface with a damp paper towel. As a final step, clean with water and wipe surface with a clean, dry paper towel.

CLEANING ALUMINIUM TOPS:

For simple dust and dirt, clean the aluminium top by using a damp cloth with soap and water. For more stubborn deposits, try using a flat edge wooden spatula to scrape off as much as possible. For more stubborn stains, try using a couple of tablespoons of white vinegar to two pints of water and mix well. Dip a clean cloth into the mixture and gently rub the exterior of the aluminium surface. Generally, it is not a good idea to use abrasive pads or cleaners on aluminium,

as the metal will scratch easily. If you must use some type of abrasive, try applying baking soda to the surface and then rubbing with a moist cloth. This will work as well as most scouring pads and is less likely to create deep scratches in the surface. Be careful not to use steel wool or scouring pads as they can leave the aluminium riddled with little scratches that make it harder to clean in the future. If you feel you must use steel wool, use the finest grade you can find and use as sparingly as possible with as little pressure as possible. Go with the grain rather than using circular motions.

INTENDED USE

These Hotplate/Stirrers/Hotplate-Stirrers are intended for general laboratory use.

EQUIPMENT DISPOSAL

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

ENVIRONMENTAL CONDITIONS

Operating Conditions: Indoor use only.

Temperature:	5 to 40°C
Humidity:	20% to 80% relative humidity, non-condensing
Altitude:	2000 M above sea level

Non-Operating Storage:

Temperature:	-20 to 65°C
Humidity:	20% to 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

SAFETY INSTRUCTIONS

Please read the entire instruction manual before operating the Hotplate/Stirrer/Hotplate-Stirrer.

WARNING! DO NOT use the Hotplate/Stirrer/Hotplate-Stirrer in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the top plate.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly after the unit has cooled down. **DO NOT** immerse the unit for cleaning. Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.

CAUTION! The top plate can reach 500°C, **DO NOT** touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.

WARNING! Units are **NOT** explosion proof. Use caution when heating volatile materials.

Replace the ceramic top immediately if damaged. A damaged top can break in use.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground - Protective Conductor Terminal

Alternating Current

STANDARDS & REGULATIONS

VWR International hereby declares under its sole responsibility that the construction of this product conforms in accordance with the following standards:

Safety standards:

- IEC 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use. Part: General Requirements.
- IEC 61010-2-010 Part II: Particular requirements for laboratory equipment for the heating of materials.
- IEC 61010-2-051 Part II: Particular requirements for laboratory equipment for mixing and stirring.
- UL Std. No. 61010-1

EMC standards:

- | | |
|-------------------|-----------------|
| EN61326-1 Class A | EN61000-3-3/3-2 |
| EN6100-4-5 | EN61000-4-4 |
| EN55022-B | EN61000-4-3 |
| EN61000-4-11 | EN61000-4-6 |

Associated EU guidelines:

- EMC directive 2004/108/EC
LVD directive 2006/95/EC

10 x 10 HOTPLATE/STIRRER/HOTPLATE-STIRRER SPECIFICATIONS

Dimensions (L x W x H): 27,4 x 16,7 x 10,8cm

Top plate dimensions (L x W): 10,2 x 10,2cm

Electrical (50/60 Hz):

Hotplate	Stirrer	Hotplate-Stirrer
1,5 amps/350 watts	0,2 amps/50 watts	1,7 amps/400 watts

Fuses: 5mm x 20mm, 5 amp quick acting

Temperature range:

ceramic: ambient +5°C to 500°C

aluminium: ambient +5°C to 400°C

Temperature stability: ceramic*: +/-3%

aluminium*: +/-2%

Speed range: 60 to 1600rpm

Speed stability: +/-2%

Capacity: 600ml, gross weight should not exceed 9kg

Controls: see diagrams

Ship weight: 2,8kg

* Below 100°C +/-2°C. Environmental and sample conditions permitting.

NOTE: On Advanced units, the max. temperature setting on the display is 500°C for a ceramic top and 400°C for an aluminium top.

Plug	Stirrer - Ceramic	Hotplate - Ceramic	Hotplate-Stirrer - Ceramic
EU	444-0566	444-0575	444-0593
UK	444-0567	444-0576	444-0594
CH	444-0568	444-0577	444-0595
Plug	Hotplate - Aluminium	Hotplate-Stirrer - Aluminium	
EU	444-0578	444-0596	
UK	444-0579	444-0597	
CH	444-0580	444-0598	

ADVANCED SERIES 10 X 10

VWR Advanced Series 10 x 10,
Ceramic Hotplate-Stirrer

caution hot top
indicator light

heat adjust knob

recall last temp. button


temperature display

stir adjust knob

ADVANCED SERIES 18 X 18

VWR Advanced Series 18 x 18,
Ceramic Hotplate-Stirrer

18 X 18 HOTPLATE/STIRRER/HOTPLATE-STIRRER SPECIFICATIONS

Dimensions (L x W x H): 37,5 x 25 x 10,8cm

Top plate dimensions (L x W): 17,8 x 17,8cm

Electrical (50/60 Hz):

	Hotplate	Stirrer	Hotplate-Stirrer
	4,4 amps/1000 watts	0,2 amps/50 watts	4,6 amps/1050 watts

Fuses:	Hotplate	Stirrer	Hotplate-Stirrer
	HP / HPS Stirrer	5mm x 20mm, 16 amp slow blow 5mm x 20mm, 5 quick acting	

Temperature range:	ceramic:	aluminium:
	ambient +5°C to 500°C	ambient +5°C to 400°C

Temperature stability:	ceramic*:	aluminium*:
	+/-3%	+/-2%

Speed range: 60 to 1600rpm

Speed stability: +/-2%

Capacity: 2500ml, gross weight should not exceed 18,1kg

Controls: see diagrams

Ship weight: 4,8kg

* Below 100°C +/-2°C. Environmental and sample conditions permitting.

NOTE: On Advanced units, the max. temperature setting on the display is 500°C for a ceramic top and 400°C for an aluminium top.

Plug	Stirrer - Ceramic	Hotplate - Ceramic	Hotplate-Stirrer - Ceramic
EU	444-0569	444-0581	444-0599
UK	444-0570	444-0582	444-0600
CH	444-0571	444-0583	444-0601
Plug	Hotplate - Aluminium	Hotplate-Stirrer - Aluminium	
EU	444-0584	444-0602	
UK	444-0585	444-0603	
CH	444-0586	444-0604	

25 x 25 HOTPLATE/STIRRER/HOTPLATE-STIRRER SPECIFICATIONS

Dimensions (L x W x H): 45,5 x 33,0 x 10,8cm

Top plate dimensions (L x W): 25,4 x 25,4cm

Electrical (50/60 Hz):

Hotplate	Stirrer	Hotplate-Stirrer
6,7 amps/1550 watts	0,2 amps/50 watts	7,0 amps/1600 watts

Fuses:	HP / HPS Stirrer	5mm x 20mm, 16 amp slow blow 5mm x 20mm, 5 quick acting
--------	---------------------	--

Temperature range:	ceramic:	ambient +5°C to 500°C
	aluminium:	ambient +5°C to 400°C

Temperature stability:	ceramic*:	+/-3%
	aluminium*:	+/-2%

Speed range: 60 to 1600rpm

Speed stability: +/-2%

Capacity: 6000ml, gross weight should not exceed 22,6kg

Controls: see diagrams


Ship weight: 7,2kg

* Below 100°C +/-2°C. Environmental and sample conditions permitting.

NOTE: On Advanced units, the max. temperature setting on the display is 500°C for a ceramic top and 400°C for an aluminium top.

Plug	Stirrer - Ceramic	Hotplate - Ceramic	Hotplate-Stirrer - Ceramic
EU	444-0572	444-0587	444-0605
UK	444-0573	444-0588	444-0606
CH	444-0574	444-0589	444-0607
Plug	Hotplate - Aluminium	Hotplate-Stirrer - Aluminium	
EU	444-0590	444-0608	
UK	444-0591	444-0609	
CH	444-0592	444-0610	

ADVANCED SERIES 25 X 25

VWR Advanced Series 25 x 25, Aluminium Hotplate-Stirrer with optional Probe Kit and glassware (glassware not included)

HEATING OPERATING INSTRUCTIONS

The Hotplates and Hotplate-Stirrers have a micro-processor controlled heater that is designed to bring samples to temperature quickly and accurately.

1. Getting ready:

- Turn the heat knob to the off position. Plug power cord into a properly grounded 3-prong outlet.
- Place a vessel with solution and the appropriate accessories in the center of the top plate. This is important because the vessel should be over the hottest part of the top plate.

2. Setting temperature for Advanced Series :

- Turn the heat knob clockwise until the display reaches the desired heat setting. The display will flash the set-point temperature until the temperature is reached, at which time the display will stopflashing. When the heat is turned on the indicator light above the heat knob is illuminated. Removing or adding more to a sample content could cause the temperature to fluctuate. If this occurs, the display will againstart to flash until the set-point value is stabilized.
- Temperature adjustments can be made without interrupting heating by turning the heat knob clockwise to increase heat or counter-clockwise to decrease heat.
- To stop heating, turn the heat knob to the off position. Your vessel can then be removed.

Caution hot top indicator:

The caution hot top indicator light warns that the top plate is too hot to touch. The caution hot top indicator light will illuminate when the heat is turned on and remain on until the top plate cools down.

Recall last temperature button:

Advanced Series Hotplate/Hotplate-Stirrers have a built-in memory that allows users to recall the last set temperature, even after the unit has been turned off. (The temperature in memory is the last temperature that ran for more than 5 minutes.)

Max Temperature:

If the unit has an aluminium top, the max temperature setting is 400°C. If the unit has a ceramic top, the max temperature setting is 500°C.

OPERATING TIPS

The unit may overshoot the temperature up to 10°C before stabilizing at the set-point. The two methods suggested to minimize overshoot are:

- Metal containers minimize overshoot. **CAUTION:** When heating metal containers on a ceramic top it is recommended to use the lowest temperature setting possible to limit thermal stress to the ceramic top.
- If a glass vessel is used, anticipate overshoot. Start with a setting 5-10°C below the desired temperature. When the temperature stabilizes at this lower setting, turn the heat knob to the final temperature. Overshoot is then reduced to about 1°C.

The temperature display on the advanced units show the actual temperature of the heater not the top plate or sample. The vessel contents being heated may be at a lower temperature depending on the size and insulating qualities of the vessel. It may be beneficial to monitor the temperature of the vessel contents and adjust the set-point temperature accordingly.

TYPICAL TIME TO BOIL WATER

The chart below is an example of an approximate time to boil for the specified amount of water in a specific vessel.

These values are only approximate and can vary from unit to unit. Values are based on 23°C water in an ambient environment of 23°C.

CERAMIC TOP			ALUMINIUM TOP		
Size of Top Plate	Volume of Water	Time	Size of Top Plate	Volume of Water	Time
10 x 10	500ml in 600ml beaker	40	10 x 10	500ml in 600ml beaker	40
18 x 18	1500ml in 2,000ml beaker	45	18 x 18	1500ml in 2000ml beaker	60
25 x 25	5000ml in 6000ml Erlenmeyer Flask	60	25 x 25	5000ml in 6000ml Erlenmeyer Flask	70

TEMPERATURE CALIBRATION PROCEDURE (SINGLE POINT CALIBRATION, ADVANCED SERIES)

This method can be used for calibrating the top plate surface of the Advanced Series unit.

To set a Single Point Calibration (SPC) adjustment:

- a. Turn the heat knob clockwise until the display reaches the desired heat setting, and let the unit stabilize to the user input temperature. The heater temperature has stabilized when the temperature display is no longer blinking and the unit will beep five times.
- b. Wait ten minutes for the surface temperature to stabilize. Measure the top plate with a traceable surface temperature measuring device.
- c. Press and hold the Last Temperature button. The display will start to scroll through the available SPC options (“UP”, “dn”, “SEt” and “dEL”) as long as the Last Temperature button is held down. Once you select “UP” or “dn” mode and release the Last Temperature button the unit will beep two times and the display will begin to blink quickly. This lets you know you are programming in SPC mode.

NOTE: There is a thirty second time out (Last Temperature button must be pressed at least once in thirty seconds, or unit will exit SPC mode). **Do not touch the heat knob while in SPC mode.** If the heat knob is adjusted during this procedure you will exit SPC mode.

- d. Once the desired option is displayed, release the Last Temperature button. Please see the explanation below for each option.
- e. Selecting the “SEt” option saves the Single Point Calibration adjustment for that temperature set point and allows you to exit this SPC mode (see Section c). When the “SEt” option is selected “SEt” will be displayed. To save the current SPC point and exit the SPC programming mode, release the Last Temperature button when “SEt” is on the display. The display will now show your set point temperature with a decimal point for that setting.
- f. Select the “UP” option if your externally measured temperature of the top plate is **higher** than the set point on the display. When the “UP” option is selected the current SPC adjusted temperature is displayed and blink-

ing quickly. **To increase the SPC, press and release the Last Temperature button multiple times (do not touch heat knob)** until the display reads the value you recorded as the measured temperature of the top plate. Changes are not saved until the “SEt” option is selected (and the Last Temperature button is pressed and released), if the temperature is adjusted too high, delete the SPC adjustment and repeat procedure.

- g. Select the “dn” option if your externally measured temperature of the top plate is **lower** than the set point on the display. When the “dn” option is selected the current SPC adjusted temperature will be displayed and blinking quickly. **To decrease the SPC, press and release the Last Temperature button multiple times (do not touch the heat knob)** until the display reads the value you recorded as the measured temperature of the top plate. Changes are not saved until the “SEt” option is selected (and the Last Temperature button is pressed and released), if the temperature is adjusted too low, delete the SPC adjustment and repeat procedure.
- h. Selecting the “dEL” option will delete **all** Single Point Calibration points and allow you to exit this SPC mode (see Section c). When the “dEL” option is selected “dEL” will be displayed. To delete **all** SPC points and exit the SPC mode release the Last Temperature button when “dEL” is on the display.
- i. For set point temperatures with a SPC adjustment, there will be a decimal point in the display. Once the SPC adjustment is set, the display will blink while the unit’s temperature is settling. When the SPC set point is reached, the display will stop blinking and the unit will beep five times.
- j. This process may be repeated for up to **three separate** set points. If a fourth SPC set point is entered, the first set point will be overwritten. To readjust an existing SPC set point, you must delete the current settings (all SPC points will be deleted, and the decimal points will no longer be displayed at those temperatures) and repeat the SPC procedure. If SPC adjustments are **not deleted** prior to resetting SPC for a set point then the temperature adjustment will not be accurate
- k. The SPC adjustments are limited to the maximum and minimum temperatures and limits allowed by the particular unit.

STIRRING OPERATING INSTRUCTIONS

The micro-processor controlled ramping feature slowly increases speed until the set-point is reached. This feature helps to avoid splashing, improves magnetic coupling and provides excellent low end control. The micro-processor also monitors and regulates the stirring speed, sensing your requirements whether you're stirring an aqueous, viscous or semi-solid solution.

Initial stirring speed may exceed set speed if the following conditions exist:

1. The stirrer is set at a low speed and the stirrer has not been operated for an extended period of time.
2. The stirrer is set at a low speed and it is the stirrer's initial use.

1. Getting ready:

- a. Turn stir knob to the off position. Plug power cord into a properly grounded outlet.
- b. Place a vessel with solution and the appropriate spin bar in the center of the top plate.

2. Setting speed:

- a. Turn the stir knob clockwise until the pointer reaches the desired speed setting. The stir indicator light above the stir knob will illuminate to indicate the stirring feature is in use. The stir indicator light will blink while reaching the set-point. Once the set-point is reached the light will remain lit.
- b. Speed adjustments can be made without interrupting stirring by turning the stir knob clockwise to increase speed, or counter-clockwise to decrease speed.
- c. To stop stirring, turn the stir knob to the off position. Your vessel can then be removed.

Stir protection for Hotplate-Stirrers:

If stirrer motor stops or fails, the unit will automatically shut down the heater.

OPERATING TIPS

The stirrer increases speed at a steady rate until the set-point is reached, if the stir bar is too large or the liquid is too viscous, the stirrer may not reach its set-point. The set-point speed needs to be reduced. The magnetic strength of stir bars reduce over time and may need to be replaced.

Stirring vessels in oil baths:

When heating and stirring a reaction vessel within an oil bath or similar set-up, the stirring function will stir up to approximately 2,54cm from the top plate. The stirring speed will vary according to liquid viscosity, stir bar length and distance from the top plate. Adjust one or all of these to achieve the desired stirring speed.

EXAMPLE: The closer the reaction vessel is to the top plate the stronger the magnetic connection.

TECHNICAL SERVICE

For information or technical assistance contact your local VWR representative or visit vwr.com.

TROUBLESHOOTING

If the unit gives an error code, immediately switch the unit off. See error table below for proper corrective action. If the error cannot be cleared, please contact your VWR representative for repairs.




Error	Cause of Error	How to Fix
E1	An "E1" error means the heater sensor is open or malfunctioned.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E2	An "E2" error means the heater sensor shorted or malfunctioned.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E3	An "E3" error means there is either no motion on the motor or the motor is not working properly.	Reset the unit by rotating the knob for speed to the off position until it clicks then turn it back on. (Also the heater will shut off. If you were heating, rotate the knob for heat to the off position then back on again if you were heating). If it still doesn't work, please contact your VWR representative for repair.
E6	An "E6" error means there is an internal electronics system error.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E8	An "E8" error means the unit had a catastrophic over temperature condition or temperature runaway condition (temperature greater than 600°C) and therefore automatically shut down to prevent damage.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E9	An "E9" error means the heater failed. This might occur if heater temperature fails to rise when asked to, or there is a sudden drop in heater temperature for no apparent reason.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.

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