Instruction Manual

Professional Series

Stirrers

Hotplates

Hotplate-Stirrers

TABLE OF CONTENTS

Package Contents	1
Warranty	1
Installation	2
Maintenance & Servicing	2
Intended Use	2
Equipment Disposal	2
Environmental Conditions	2
Safety Instructions	3
Standards & Regulations	3
Control Panel	4
18 x 18 Specifications	5
25 x 25 Specifications	6
Heating Operating Instructions	7-10
Stirring Operating Instructions	11
Function Mode	12-14
Technical Service	15
Troubleshooting	15-16
Replacement Parts	17-18

PACKAGE CONTENTS

- Hotplate, Stirrer or Hotplate-Stirrer
- 234cm detachable power cord
- Spin bar (stirring models only)
- Probe Kit (heating models only): 20.3cm stainless steel probe, 45.7cm support rod, thermometer/temperature probe extension clamp, hook connector
- Instruction manual

WARRANTY

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of purchase. If a defect is present, VWR will, at its option, repair, replace, or refund the purchase price of this product at no charge to you, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. For your protection, items being returned must be insured against possible damage or loss. This warranty shall be limited to the replacement of defective products. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

INSTALLATION

Upon receiving the VWR Hotplate/Stirrer/Hotplate-Stirrer, check to ensure that no damage has occurred in shipment. It is important that any damage that occurred in transport is detected at the time of unpacking. If you do find such damage the carrier must be notified immediately.

After unpacking, place the Hotplate/Stirrer/Hotplate-Stirrer on a level bench or table, away from explosive vapors. Ensure that the surface on which the unit is placed will withstand typical heat produced by the unit and place the unit a minimum of 15,2 centimeters from vertical surfaces. Always place the unit on a sturdy work surface.

The Hotplate/Stirrer/Hotplate-Stirrer is supplied with a power cord that is inserted into the IEC connector on the back of the unit first, then it can be plugged into a properly grounded outlet. The unit plugs into a 230 volt, 50/60 Hz source.

MAINTENANCE & SERVICING

The Hotplate/Stirrer/Hotplate-Stirrer is built for long, trouble-free, dependable service. It needs no user maintenance beyond keeping the surfaces clean. The unit should be given the care normally required for any electrical appliance. Avoid wetting or unnecessary exposure to fumes. Spills should be removed promptly after the unit has cooled down. Do not use a cleaning agent or solvent on the front panel or top plate which is abrasive or harmful to plastics, nor one which is flammable. Always ensure the power is disconnected from the unit prior to any cleaning. If the unit ever requires service, contact your VWR representative.

CLEANING CERAMIC TOPS:

First remove any burnt-on deposits or spills from the top plate with a scraper (similar to scraping paint off of windowpanes in your home). For your safety, please wear an insulated mat when using a metal scraper. When the top plate has cooled, apply a few dabs of a non-abrasive cleaner over the surface with a damp paper towel. As a final step, clean with water and wipe surface with a clean, dry paper towel.

CLEANING ALUMINIUM TOPS:

For simple dust and dirt, clean the aluminium top by using a damp cloth with soap and water. For more stubborn deposits, try using a flat edge wooden spatula to scrape off as much as possible. For more stubborn stains, try using a couple of tablespoons of white vinegar to two pints of water and mix well. Dip a clean cloth into the mixture and gently rub the exterior of the aluminium

surface. Generally, it is not a good idea to use abrasive pads or cleaners on aluminium, as the metal will scratch easily. If you must use some type of abrasive, try applying baking soda to the surface and then rubbing with a moist cloth. This will work as well as most scouring pads and is less likely to create deep scratches in the surface. Be careful not to use steel wool or scouring pads as they can leave the aluminium riddled with little scratches that make it harder to clean in the future. If you feel you must use steel wool, use the finest grade you can find and use as sparingly as possible with as little pressure as possible. Go with the grain rather than using circular motions.

INTENDED USE

These Hotplate/Stirrers/Hotplate-Stirrers are intended for general laboratory use.

EQUIPMENT DISPOSAL

This equipment must not be disposed of with unsorted waste. It is your responsibility to correctly dispose of the equipment at life-cycle-end by handing it over to an authorized facility for separate collection and recycling. It is also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect the persons involved in the disposal and recycling of the equipment from health hazards.

For more information about where you can drop off your waste of equipment, please contact your local dealer from whom you originally purchased this equipment. By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

ENVIRONMENTAL CONDITIONS

Operating Conditions: Indoor use only.

Temperature:	5 to 40°C
Humidity:	20% to 80% relative humidity, non-condensing
Altitude:	2000 M above sea level

Non-Operating Storage:

Temperature:	-20 to 65°C
Humidity:	20% to 80% relative humidity, non-condensing

Installation Category II and Pollution Degree 2 in accordance with IEC 664.

SAFETY INSTRUCTIONS

Please read the entire instruction manual before operating the Hotplate, Stirrer, Hotplate-Stirrer.

WARNING! DO NOT use the Hotplate, Stirrer, Hotplate-Stirrer in a hazardous atmosphere or with hazardous materials for which the unit was not designed. Also, the user should be aware that the protection provided by the equipment may be impaired if used with accessories not provided or recommended by the manufacturer, or used in a manner not specified by the manufacturer.

Always operate unit on a level surface for best performance and maximum safety.

DO NOT lift unit by the top plate.

CAUTION! To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.

Spills should be removed promptly after the unit has cooled down. **DO NOT** immerse the unit for cleaning. Alkalis spills, hydrofluoric acid or phosphoric acid spills may damage the unit and lead to thermal failure.

CAUTION! The top plate can reach 500°C, **DO NOT** touch the heated surface. Use caution at all times. Keep the unit away from explosive vapors and clear of papers, drapery and other flammable materials. Keep the power cord away from the heater plate.

DO NOT operate the unit at high temperatures without a vessel/sample on the top plate.

WARNING! Units are **NOT** explosion proof. Use caution when heating volatile materials.

Replace the ceramic top immediately if damaged. A damaged top can break in use.

DO NOT operate the unit if it shows signs of electrical or mechanical damage.

Earth Ground - Protective Conductor Terminal

Alternating Current

STANDARDS & REGULATIONS

VWR International hereby declares under its sole responsibility that the construction of this product conforms in accordance with the following standards:

Safety standards:


IEC 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use. Part: General Requirements.
IEC 61010-2-010	Part II: Particular requirements for laboratory equipment for the heating of materials.
IEC 61010-2-051	Part II: Particular requirements for laboratory equipment for mixing and stirring.
UL Std. No. 61010-1	

EMC standards:

EN61326-1 Class A	EN61000-3-3/3-2
EN6100-4-5	EN61000-4-4
EN55022-B	EN61000-4-3
EN61000-4-11	EN61000-4-6

Associated EU guidelines:

EMC directive 2004/108/EC
LVD directive 2006/95/EC

CONTROL PANEL

The front panel of the Professional Series Hotplate/Stirrer/Hotplate-Stirrer contains all the controls and displays needed to operate the unit.

A. Standby button/standby indicator light: The standby indicator light will illuminate when the unit is plugged in. The unit will be in standby mode. Press the standby button to illuminate the temperature, speed, and time displays.

B. Temperature display: Temperature display: Displays the actual/setpoint temperatures in conjunction with the actual/setpoint indicator lights. **C.** Up/down arrows for set-point control. On/off button starts/stops the heating function.

D. Heat indicator light: Illuminates when the external temperature probe is not being used. The actual temperature displayed is the heater temperature.

E. Actual indicator light: Illuminates when the temperature displayed is the actual temperature of the heater/temperature probe.

F. Probe indicator light: Illuminates when the external temperature probe is plugged in. The actual temperature displayed is the probe temperature.

G. Set-point indicator light: Illuminates when the set-point temperature is displayed.


H. Caution hot top indicator light: Illuminates when the heat is turned on and warns the top plate is too hot to touch. Remains on until the top plate has cooled down.

I. Speed display: Displays the speed of the unit. **J.** Up/down arrows for set-point control. On/off button starts/stops the stirring function.

K. Time display: Displays accumulated time (continuous mode) or how much time is remaining (timed mode). The display range is from 0 to 9,999 minutes in one (1) second increments. The display will indicate minutes and seconds until the timer reaches 99 minutes and 59 seconds (99:59), then the display will automatically display minutes up to 9999. **L.** Up/down arrows for set-point control.

PROFESSIONAL SERIES 18 X 18

VWR Professional Series
18 x 18, ceramic
Hotplate-Stirrer with
Probe Kit and
glassware*

*glassware not included

18 x 18 HOTPLATE/STIRRER/HOTPLATE-STIRRER SPECIFICATIONS

Dimensions (L x W x H): 37,5 x 25 x 10,8cm

Top plate dimensions (L x W): 17,8 x 17,8cm

Electrical (50/60 Hz):

Hotplate	Stirrer	Hotplate-Stirrer
4,4 amps/1000 watts	0,2 amps/50 watts	4,6 amps/1050 watts

Fuses: 5mm x 20mm, 10 amp quick acting

Temperature range:
aluminium: ambient +5°C to 400°C
ceramic: ambient +5°C to 500°C

Temperature stability: +/-1%

Speed range: 60 to 1600rpm

Speed stability: +/-2%

Capacity: 2500ml, gross weight should not exceed 18,1kg

Controls: see page 4

Ship weight: 2,5kg

* Below 100°C +/-1°C. Environmental and sample conditions permitting.

Plug	Stirrer - Ceramic	Hotplate - Ceramic	Hotplate-Stirrer - Ceramic
EU	444-0611	444-0617	444-0629
UK	444-0612	444-0618	444-0630
CH	444-0613	444-0619	444-0631
Plug	Hotplate - Aluminium		Hotplate-Stirrer - Aluminium
EU	444-0620		444-0632
UK	444-0621		444-0633
CH	444-0622		444-0634

25 x 25 HOTPLATE/STIRRER/HOTPLATE-STIRRER SPECIFICATIONS

Dimensions (L x W x H): 45,5 x 33,0 x 10,8cm

Top plate dimensions (L x W): 25,4 x 25,4cm

Electrical (50/60 Hz):

Hotplate	Stirrer	Hotplate-Stirrer
6,7 amps/1550 watts	0,2 amps/50 watts	7,0 amps/1600 watts

Fuses: 5mm x 20mm, 16 amp quick acting

Temperature range: **aluminium:** ambient +5°C to 400°C
ceramic: ambient +5°C to 500°C

Temperature stability: +/-1%

Speed range: 60 to 1600rpm

Speed stability: +/-2%

Capacity: 6000ml, gross weight should not exceed 22,7kg

Controls: see page 4

Ship weight: 3,6kg

** Below 100°C +/-1°C. Environmental and sample conditions permitting.*

Plug	Stirrer - Ceramic	Hotplate - Ceramic	Hotplate-Stirrer - Ceramic
EU	444-0614	444-0623	444-0635
UK	444-0615	444-0624	444-0636
CH	444-0616	444-0625	444-0637
Plug		Hotplate - Aluminium	Hotplate-Stirrer - Aluminium
EU		444-0626	444-0638
UK		444-0627	444-0639
CH		444-0628	444-0640

PROFESSIONAL SERIES 25 X 25

VWR Professional Series
25 x 25, aluminium
Hotplate-Stirrer

HEATING OPERATING INSTRUCTIONS

The Hotplates and Hotplate-Stirrers have a micro-processor controlled heater that is designed to bring samples to temperature quickly and accurately.

1. Getting ready:

- Plug power cord into a properly grounded 3-prong outlet.
- Place a vessel with solution and the appropriate accessories in the center of the top plate. This is important because the vessel should be over the hottest part of the top plate.
- Press the standby button to move the unit from standby mode. The temperature and time displays will illuminate. On power-up, allow the unit approximately 30 seconds to stabilize before accurate temperature readings are displayed.

2. Setting temperature:

- Press the up/down arrows below the temperature display until you reach the desired temperature. Press the on/off button below the temperature display to start the heating function. The indicator light below the temperature display will illuminate to indicate the heating function is in use. Audible beeps will sound to indicate the set-point temperature has been reached. (LED display will toggle between actual and set point temperatures. Indicator lights will toggle between actual and set point to indicate what is displayed).
- Set-point temperature adjustments can be made without interrupting heating using the up/down arrows below the temperature display. Removing or adding more to a sample content could cause the temperature to fluctuate.
- To stop heating function, press the on/off button below the temperature display. Your vessel can then be removed.

3. Setting temperature with probe:

- Probe must be connected to the unit before turning on the heating function.
- Follow step 2a through 2c.

Caution hot top indicator:

The caution hot top indicator light warns that the top plate is too hot to touch. The caution hot top indicator light will illuminate when the heat is turned on and remain on until the top plate cools down.

4. Set timed mode:

- Press the up/down arrows below the time display until you reach the desired time.
- When the unit is started, by pressing the on/off button below the temperature or speed display, the unit will run for the selected time. There are audible beeps to indicate the time down function is complete. The unit will shut off automatically.
- To repeat for the same time, simply depress the on/off button again.
- To interrupt an automatic timing cycle before it is completed, press the on/off button below the temperature or speed display.

5. Set time to continuous mode:

- Set timer to zero (0) using the up/down arrows below time display.
- When the unit is started, by pressing the on/off button below the temperature or speed display, the time display will show actual running time.
- To stop, press the on/off button below the temperature or speed display.

Max Temperature:

If the unit has an aluminium top, the max temperature setting is 400°C. If the unit has a ceramic top, the max temperature setting is 500°C. The only exception is if the unit is a Professional Series with a ceramic top. In that case when using probe control, the software limits the Professional ceramic top unit to 400°C. If the probe is not plugged in, the unit can be set to 500°C.

*** CAUTION: If you are using a probe accessory with a Professional unit, make sure you don't exceed a sample temperature of 250°C, otherwise the probe tip might become damaged!**

HEATING OPERATING TIPS

Overshoot:

The unit may overshoot the temperature up to 10°C before stabilizing at the set-point. The two methods suggested to minimize overshoot are:

1. Metal containers minimize overshoot. **CAUTION!** When heating metal containers on a ceramic top it is recommended to use the lowest temperature setting possible to limit thermal stress to the ceramic top.
2. If a glass vessel is used, anticipate overshoot. Start with a setting 5 to 10°C below the desired temperature. When the temperature stabilizes at this lower setting, press the up/down arrows to the final temperature. Overshoot is then reduced to about 1°C.

The temperature display on the units show the actual temperature of the heater not the top plate or sample. The vessel contents being heated may be at a lower temperature depending on the size and insulating qualities of the vessel. It may be beneficial to monitor the temperature of the vessel contents and adjust the set-point temperature accordingly. If you need precise control a probe will be necessary.

TYPICAL TIME TO BOIL WATER

The chart below is an example of an approximate time to boil for the specified amount of water in a specific vessel.

These values are only approximate and can vary from unit to unit. Values are based on 23°C water in an ambient environment of 23°C.

CERAMIC TOP			ALUMINIUM TOP		
Size of Top Plate	Volume of Water	Time	Size of Top Plate	Volume of Water	Time
18 x 18	1500ml in 2000ml beaker	45	18 x 18	1500ml in 2000ml beaker	60
25 x 25	5000ml in 6000ml Erlenmeyer Flask	60	25 x 25	5000ml in 6000ml Erlenmeyer Flask	70

Probe Controlled Plate Temperature Limit (Flash point Setting):

(for Hotplates and Hotplate-Stirrers)

See Function Mode section for instructions for the Plate Limit Setting procedure for limiting the top plate temperature when not using a probe.

To adjust the maximum probe controlled plate temperature

1. Place unit in standby mode and ensure probe is properly connected.
2. Press and hold the temperature on/off and time down buttons.
3. While pressing and holding the temperature on/off and time down buttons, press and release the standby button.
 - a. The probe LED will be on during this step.
 - b. If maximum probe controlled plate temperature (flash point) is set below the maximum, the temperature LED will be off for the first 1 second of the actual and set temperature display while heating.
 - c. Use the temperature up and temperature down buttons to adjust the maximum probe controlled plate temperature (flash point).
 - d. The minimum value for this temperature is 60°C for both aluminium top and ceramic top units.
 - e. The maximum value for this temperature is 400°C for aluminium top units and 430°C for ceramic top units.
 - f. Press and release the standby button to return to standby mode.

NOTE: The indicator light below the temperature display will blink every 3 to 4 seconds to indicate this function is enabled.

EXAMPLE:

Change the plate over-temperature limit from the current 430°C to 200°C.

1. Start with the unit turned off.
2. Press and hold the on/off button below the temperature display and the time down arrow while turning the unit on, by pressing the standby button. The temperature displayed is the current plate over-temperature limit setting, 430°C.
3. Press the up/down arrows until the display reads 200°C. After approximately 2 seconds, the display will blink off, then on again to store the new setting.
4. The new over-temperature limit is now set for 200°C. Turn the unit off before using.

TEMPERATURE CALIBRATION PROCEDURE (SINGLE POINT CALIBRATION)

This method can be used for calibrating the RTD probe if connected to the unit or the top plate surface of the unit if the RTD probe is not connected.

1. If unit is in standby mode press the standby button to take unit out of standby mode.
2. To set a Single Point Calibration adjustment:
 - a. Set the unit to the desired temperature, turn the temperature control system on and let the system stabilize to the user input temperature.
 - b. While pressing and holding the standby button, press and release the temperature up button.
 - c. The unit will beep two times (if not in silent mode) and the temperature display will blink to indicate a Single Point Calibration adjustment is in progress.
 - d. Measure the sample with a traceable thermometer or temperature probe.
 - e. Use the temperature up and down buttons to adjust the display to read the actual measured temperature indicated on the external temperature measuring device. (Example: Desired temperature is 80°C. Set unit at 80°C per step (a). Follow steps (b & c). Display reads 80°C and the external temperature measuring device reads higher at 84°C. Push the temperature up arrow so the display will match the external temperature device and also read 84°C. By doing this a biased offset for 80°C will be used any time the unit is set to 80°C)
 - f. Press the standby button to save the Single Point Calibration adjustment. The temperature display will stop blinking and the decimal point will flash periodically to indicate that a Single Point Calibration adjustment is in use. This indication will occur whenever the temperature is adjusted to the temperature used in part “a” of this procedure.

- g. This process may be repeated for up to three separate set points. If a fourth calibration set point is entered, the first set point will be overwritten.
- h. The adjustments are limited to the maximum and minimum temperature range of the particular model.

The unit will now use the biased offset for that specific temperature setting and increase or decrease temperature accordingly to bring the unit temperature to set-point. The decimal point of the display will flash to indicate a biased offset is being used. All other temperature settings will use the factory default settings. The offset will be stored in memory and retained until reset.

3. To clear all existing Single Point Calibration adjustments:
 - a. While pressing and holding the standby button, press and release the temperature down button.
 - b. The unit will beep two times (if not in silent mode) to indicate all Single Point Calibration adjustments have been cleared.

STIRRING OPERATING INSTRUCTIONS

The micro-processor controlled ramping feature slowly increases speed until the set-point is reached. This feature helps to avoid splashing, improves magnetic coupling and provides excellent low end control. The micro-processor also monitors and regulates the stirring speed, sensing your requirements whether you're stirring an aqueous, viscous or semi-solid solution.

Initial stirring speed may exceed set speed if the following conditions exist:

1. The stirrer is set at a low speed and the stirrer has not been operated for an extended period of time.
2. The stirrer is set at a low speed and it is the stirrer's initial use.

1. Getting ready:

- a. Plug power cord into a properly grounded outlet.
- b. Place a vessel with solution and the appropriate accessories in the center of the top plate. This is important because the spin bar will be well within the magnetic field.
- c. Press the standby button to move the unit from standby mode. The speed and time displays will illuminate.

2. Setting speed:

- a. Press the up/down arrows below the speed display until you reach the desired speed.
- b. To start stirring function, press the on/off button below the speed display.
- c. The indicator light below the speed display will illuminate to indicate the stirring function is in use.
- d. Speed adjustments can be made without interrupting stirring by using the up/down arrows below the speed display.
- e. To stop stirring function, press the on/off button below the speed display.

3. Set timed mode:

- a. Press the up/down arrows below the time display until you reach the desired time.
- b. When the unit is started, by pressing the on/off button below the temperature or speed display, the unit will run for the selected time. There are audible beeps to indicate the time down function is complete. The unit will shut off automatically. To repeat for the same time, simply depress the on/off button again.

- c. To interrupt an automatic timing cycle before it is completed, press the on/off button below the temperature or speed display.

4. Set time to continuous mode:

- a. Set timer to zero (0) using the up/down arrows below time display.
- b. When the unit is started, by pressing the on/off button below the temperature or speed display, the time display will show actual running time.
- c. To stop, press the on/off button below the temperature or speed display.

Stir protection for Hotplate-Stirrers:

If stirrer motor stops or fails, the unit will automatically shut down the heater.

OPERATING TIPS

The stirrer increases speed at a steady rate until the set-point is reached, if the spin bar is too large or the liquid is too viscous, the stirrer may not reach its set-point. The set-point speed needs to be reduced. The magnetic strength of stir bars reduce over time and may need to be replaced.

Stirring vessels in oil baths:

When heating and stirring a reaction vessel within an oil bath or similar set-up, the stirring function will stir up to approximately one 2,54cm from the top plate. The stirring speed will vary according to liquid viscosity, spin bar length and distance from the top plate. Adjust one or all of these to achieve the desired stirring speed.

EXAMPLE: The closer the reaction vessel is to the top plate the stronger the magnetic connection.

FUNCTION MODE

Hotplate-Stirrers and Stirrers (Hotplate only models see page 13):

The steps listed below allow you to customize the functionality of your Hotplate-Stirrer or Stirrer unit.

Enable/Disable Function Settings

To enter into the Function Modes, perform the following steps:

1. Place unit in standby mode.
2. Press and hold the time down button.
3. While pressing and holding the time down button, press and release the standby button.

The speed up/down arrows allow you to toggle between functions. The time up/down arrows allow you to turn the function off or on. Press the standby button at any time to exit the Function Mode and save the current settings.

BEEP DISABLE (this function will disable the beep, except one caused by an error). Speed display will show "bEEP" and time display will show current setting "on" or "off".

- a. Press the time up or down button to toggle beeper on/off.
- b. Press the speed up button to select the next option or speed down button to select the previous option.

PLATE LIMIT SETTING (this procedure will limit the top plate temperature when the probe is **NOT** being used). Speed display will show "SEtL" and time display will show current setting (degrees C).

- a. This option will be skipped on Professional Stirrer models.
- b. The heat LED will be on during this step.
- c. Press the time up or down button to adjust the current plate temperature limit.
- d. Press the speed up button to select the next option or speed down button to select the previous option.

E3 (motor error) DISABLE (Disabling this function will enable a Hotplate-Stirrer to continue to heat in the event of a motor error or E3 error which could be caused by highly viscous liquids). Speed display will show "S E" and time display will show current setting "on" or "off".

- a. Press the time up or down button to toggle Speed Error (motor error) on/off.
- b. Press the speed up button to select the next option or speed down button to select the previous option.

E7 (Probe Out of Sample) DISABLE (Disabling this function will enable a unit to continue to heat a sample in the event the sample can not reach the set point. This may be due to the fact the set point is above what the sample can achieve. Speed display will show "POOS" and time display will show current setting "on" or "off").

- a. This option will be skipped on Professional Stirrer models
- b. Press the time up or down button to toggle Probe Out of Sample on/off.
- c. Press the speed up button to select the next option or speed down button to select the previous option.

POWER FAILURE AUTO-RECOVERY DISABLE

Speed display will show "PF" and time display will show current setting "on" or "off".

- a. Press the time up or down button to toggle Power Fail option on/off. On means that the unit will automatically recover from a power failure, OFF means that the unit will indicate that a power failure has occurred, but always return the unit to standby mode in the event of a power failure.
- b. Press the Speed Up button to select the next option or Speed Down button to select the previous option.

FUNCTION MODE

During normal operation (factory default setting) if power failed while heating, stirring, or timing the unit will recover and resume the operation being performed. If the unit was running in continuous mode, it will resume in continuous mode, if the unit was running in timed mode, the time will reset to the original set time and begin counting down again. Turning this feature off will disable this function. In the event of a power failure if the function is turned off, the unit will not recover, but indicate a power failure occurred by a "PF" in the time display. To clear the "PF" indication (and place the unit in standby) press and release the standby button.

TIMER DISABLE FUNCTIONS

The factory default setting allows for both heating and stirring functions to turn off simultaneously when the timer reaches zero. The Heater/Timer and Speed/Timer Disable functions below allow you to change the default setting and program what is turned off when the timer reaches zero.

HEATER/TIMER DISABLE

Speed display will show "tHEt" and time display will show current setting ("on" or "off").

- a. This option will be skipped on Professional Stirrer models
- b. Press the time up or down button to toggle Heater/Timer interaction option on/off. On means that the unit will shut off the heater when the timer counts down to zero, OFF means that the heater will continue to run when the timer counts down to zero.
- c. Press the speed up button to select the next option or speed down button to select the previous option.

SPEED/TIMER DISABLE

Speed display will show "tSPd" and time display will show current setting ("on" or "off").

- a. Press the time up or down button to toggle Speed/Timer interaction option on/off. On means that the unit will shut off the motor when the

timer counts down to zero, OFF means that the motor will continue to run when the timer counts down to zero.

- b. Press the speed up button to select the next option or speed down button to select the previous option.

Press the standby button at any time to exit the Function Mode and save current settings. Power Fail setting, Temperature Limit, and Timer functions are saved in the memory. All other options are not stored in the memory and restored to the default setting (enabled) when power is cycled.

Hotplates (Hotplate-Stirrers and Stirrers see page 12):

The steps listed below allow you to customize the functionality of your Hotplate only unit.

Enable/Disable Function Settings

To enter into the Function Modes, perform the following steps:

1. Place unit in standby mode.
2. Press and hold the time down button.
3. While pressing and holding the time down button, press and release the standby button.

The temperature up/down arrows allow you to toggle between functions. The time up/down arrows allow you to turn the function off or on. Press the standby button at any time to exit the Function Mode and save the current settings.

FUNCTION MODE

BEEP DISABLE

Temperature display will show “bEEP” and time display will show current setting (“on” or “off”).

- a. Press the time up or down button to toggle beeper on/off.
- b. Press the temperature up button to select the next option or temperature down button to select the previous option.

PLATE LIMIT SETTING (this procedure will limit the top plate temperature when the probe is NOT being used). Temperature display will show “SEtL” and time display will show current setting (degrees C).

- a. The heat LED will be on during this step.
- b. Press the time up or down button to adjust the current plate temperature limit.
- c. Press the temperature up button to select the next option or temperature down button to select the previous option.

E7 (probe out of sample) DISABLE (Disabling this function will enable a unit to continue to heat a sample in the event the sample can not reach the set point. This may be due to the fact the set point is above what the sample can achieve.). Temperature display will show “POOS” and time display will show current setting (“on” or “off”).

- a. Press the time up or down button to toggle Probe Out of Sample on/off.
- b. Press the temperature up button to select the next option or temperature down button to select the previous option.

POWER FAILURE AUTO RECOVERY DISABLE

Temperature display will show “PF” and time display will show current setting (“on” or “off”).

- a. Press the time up or down button to toggle Power Fail option on/off.
On means that the unit will automatically recover from a power failure, OFF means that the unit will indicate that a power failure has occurred, but always return to standby mode in the event of a power failure.
- b. Press the temperature up button to select the next option or temperature down button to select the previous option.

During normal operation (factory default setting) if power failed while heating, stirring, or timing the unit will recover and resume the operation being performed. If the unit was running in continuous mode, it will resume in continuous mode, if the unit was running in timed mode, the time will reset to the original set time and begin counting down again. Turning this feature off will disable this function. In the event of a power failure if the function is turned off, the unit will not recover, but indicate a power failure occurred by a “PF” in the time display. To clear the “PF” indication (and place the unit in standby) press and release the standby button.

Press the standby button at any time to exit the Function Mode and save current settings. Power Fail setting and Temperature Limit are saved in the memory. All other options are not stored in the memory and restored to the default setting (enabled) when power is cycled.

TECHNICAL SERVICE

For information or technical assistance contact your local VWR representative or visit vwr.com.

TROUBLESHOOTING

If the unit gives an error code, immediately switch the unit off. See error table below for proper corrective action. If the error cannot be cleared, please contact your VWR representative for repairs.


Error	Cause of Error	How to Fix
E1	An "E1" error means the heater sensor is open or malfunctioned.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E2	An "E2" error means the heater sensor shorted or malfunctioned.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E3	An "E3" error means there is either no motion on the motor or the motor is not working properly.	Reset the unit by pressing the on/off button for speed. (Also the heater will shut off. If you were heating, press the On/Off button for heat to reset). If it still doesn't work, please contact your VWR representative for repair.
E4	An "E4" error means there is an open probe. It is detected when the heating is started in probe control (with probe connected); after heating begins, the probe fails or is unplugged during the heating process. (If you want to run unit in PLATE control, heating must be started without external probe plugged in.)	Reset the unit by pressing the heat on/off button twice. Also make sure the probe is plugged in to the back of the unit. If it still does not work, please contact your VWR representative for repair.
E5	An "E5" error means the probe shorted or malfunctioned.	Reset the unit by pressing the heat on/off button twice. Also try plugging in another probe that is known to be working. If this doesn't work, please contact your VWR representative for repair.

Error codes continued on next page.

Error	Cause of Error	How to Fix
E6	An "E6" error means there is an internal electronics system error.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E7*	An "E7" error means the probe is not able to reach set point. It could mean the probe is out of the sample, or the temperature of the probe doesn't rise after a certain amount of time, or the probe temperature falls when the heater temperature is rising.	Reset the unit by pressing the heat on/off button twice. Also make sure probe is in the sample. If it still does not work, please contact your VWR representative for repair.
E8	An "E8" error means the unit had a catastrophic over temperature condition or temperature runaway condition (temperature greater than 600°C) and therefore automatically shuts down to prevent damage.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
E9	An "E9" error means the heater failed. This might occur if heater temperature fails to rise when asked to, or there is a sudden drop in heater temperature for no apparent reason.	This error cannot be fixed by the end user. Please contact your VWR representative for repair.
PF	AC power failure	Unit will go into standby and all functions will be shut off. Press the standby button to clear the PF error code. Then press the standby button to turn the unit back on. The heating and stirring functions must be restarted.


* When using probe control, the temperature setpoint should be adjusted to the desired sample temperature. If the temperature setpoint is higher than the sample can achieve, an E7 error will occur. Reduce sample volume or temperature setpoint value. For Example: Water has a theoretical temperature limit of 100°C (boiling). A temperature setpoint greater than 100°C will cause an E7 error.

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