

Datasheet

Pdcd11 polyclonal antibody

Catalog Number: PAB15709

Regulation Status: For research use only (RUO)

Product Description: Rabbit polyclonal antibody raised against partial recombinant Pdcd11.

Immunogen: Recombinant GST fusion protein corresponding to 170 mouse Pdcd11.

Sequence:

TKSEKYKEAGELYNRMLKRFRQEKA V WIKYGAFVLGR
SQAGASHRVLQRALECLPAKEHVDVIVKFAQLEFQLG
DVERAKAIFENTLSTYPKRTDVWSVYIDMTIKHGSQTA
VRDIFERVIIHLSLAPKRMKFFFKRYLDYEKQHGTEKDV
QAVKAKALEYVEAKSSALED

Host: Rabbit

Reactivity: Human

Applications: WB

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Specificity: Specific to recombinant protein GX0685.

This antibody detects mPDCD11 protein. It also recognizes human PDCD11 protein.

Form: Liquid

Recommend Usage: Western Blot (1:1000)

The optimal working dilution should be determined by the end user.

Storage Buffer: In PBS (50% glycerol, 0.02% sodium azide)

Storage Instruction: Store at -20°C.

Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 18572

Gene Symbol: Pdcd11

Gene Alias: 1110021I22Rik, ALG-4, Pdcd7, mKIAA0185

References:

1. Prediction of the coding sequences of mouse homologues of KIAA gene: III. the complete nucleotide sequences of 500 mouse KIAA-homologous cDNAs identified by screening of terminal sequences of cDNA clones randomly sampled from size-fractionated libraries. Okazaki N, Kikuno R, Ohara R, Inamoto S, Koseki H, Hiraoka S, Saga Y, Nagase T, Ohara O, Koga H. DNA Res. 2003 Aug 31;10(4):167-80.