

## Datasheet

### SRB4 polyclonal antibody

**Catalog Number:** PAB1164

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against full length recombinant *Saccharomyces cerevisiae* SRB4.

**Immunogen:** Recombinant protein corresponding to full length *Saccharomyces cerevisiae* SRB4.

**Host:** Rabbit

**Applications:** ELISA, WB-Ce

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Form:** Liquid

**Recommend Usage:** The optimal working dilution should be determined by the end user.

**Storage Buffer:** In buffer containing 0.1% sodium azide

**Storage Instruction:** Store at 4°C.

**Entrez GeneID:** 856743

**Gene Symbol:** SRB4

**Gene Alias:** MED17