

Product Specification Sheet

VWR CYLINDRICAL TEST TUBES WITHOUT RIM

VWR International certifies that the products comply in every respect with the relevant particular specifications and/or drawings and/or relevant standards and regulations in force, mentioned in this document.

VWR International certifie que les produits sont conformes à tous égards aux spécifications particulières appropriées et / ou dessins, normes et règlements en vigueur, mentionnées dans le présent document

VWR International bestätigt, dass das Produkte allen zutreffenden Spezifikations- und/oder Zeichnungs- und/oder Standardanforderungen sowie den im Dokument aufgeführten gültigen Vorschriften entspricht.

Cat. No.	Length (mm)	Ø ext.	Capacity	Material	Type	Pk
216-1044	100 mm	16 mm	10 ml	PS	Non graduated High thickness	2.000
216-1046	100 mm	16 mm	10 ml	PP	Graduated at 2.5 - 5 - 10ml	2.000
216-1163	55 mm	11.5 mm	3 ml	PS	Non graduated	6.000
216-1164	55 mm	11.5 mm	3 ml	PP	Non graduated	6.000
216-1165	75 mm	12 mm	3.5 ml	PP	Sorwall type bags	4.000
216-1167	75 mm	13 mm	5 ml	PP	Non graduated	4.000

Other Characteristics:

- Universal tubes with round bottom for a wide range of applications.
- Do not expose or leave near open flames, do not leave the product around as it is not biodegradable.
- IVD, CE marked.

Disclaimer: VWR wants to state that these specifications can only be regarded as general information and will not discharge the user from his obligation to ensure the product is suitable for the intended use - If there is any doubt ask us.

Date of issue: **September 2016**