




## B7-H4 Antibody [H74]

CATALOG NUMBER: 51-351



Immunohistochemistry of human pancreas tissue stained using B7-H4 Monoclonal Antibody.

### Specifications

<b>SPECIES REACTIVITY:</b>	Human
<b>TESTED APPLICATIONS:</b>	FACS, IHC
<b>APPLICATIONS:</b>	B7-H4 antibody can be used in Western Blot, immunohistochemistry starting at 10 ug/mL, and flow cytometry starting at 1:10.
<b>USER NOTE:</b>	Optimal dilutions for each application to be determined by the researcher.
<b>IMMUNOGEN:</b>	Human VTCN1
<b>HOST SPECIES:</b>	Mouse

### Properties

<b>PURIFICATION:</b>	Affinity Purified
<b>PHYSICAL STATE:</b>	Liquid
<b>BUFFER:</b>	PBS, pH 7.2, 150 mM sodium chloride, 0.09% sodium azide
<b>STORAGE CONDITIONS:</b>	Store B7-H4 antibody at 4 °C. As with all antibodies avoid freeze/thaw cycles
<b>CLONALITY:</b>	Monoclonal
<b>ISOTYPE:</b>	IgG1,k
<b>CONJUGATE:</b>	Unconjugated

### Additional Info

<b>ALTERNATE NAMES:</b>	VTCN1, B7 superfamily member 1, B7 Superfamily, member 1, B7h.5, B7S1, B7H4, PRO1291, VCTN1, RP11-229A19.4, B7 family member, H4, B7-H4, B7X, Protein B7S1
<b>ACCESSION NO.:</b>	Q7Z7D3
<b>PROTEIN GI NO.:</b>	74759262
<b>OFFICIAL SYMBOL:</b>	VTCN1
<b>GENE ID:</b>	79679

### Background

